

# bkFusion Kits

bkFusion requires UA2601 Prostate or UA2604 Urology Advanced Bundle

## TRANSRECTAL FUSION

Option	Systems	Transducer	Accessories
UA2379-01 bkFusion TR	bk3000/bk5000	E14C4t/ E13C2	<ul style="list-style-type: none"> <li>UA2372-01 bkFusion TR software</li> <li>UA2374 bkFusion license</li> <li>UA2371 EM sensor</li> <li>EM control unit and shelf</li> <li>EM transmitter &amp; stand</li> <li>Universal clamp</li> <li>Cables</li> </ul>

Note that sensor clamps for transducers must be ordered separately: UA2377 sensor clamp for E14C4t / UA2399 sensor clamp for E13C2

## TRANSPERINEAL FUSION

Option	Systems	Transducer	Accessories
UA2390-01 bkFusion TP Stabilized with stepper	bk3000/bk5000	E14CL4b	<ul style="list-style-type: none"> <li>UA2373-01 bkFusion TP software</li> <li>UA2374 bkFusion license</li> <li>UA2024 CIVCO EX3 Stepper Kit (for E14CL4b)</li> </ul>
UA2390-04 bkFusion TP Stabilized without stepper	bk3000/bk5000	E14CL4b	<ul style="list-style-type: none"> <li>UA2373-01 bkFusion TP software</li> <li>UA2374 bkFusion license</li> </ul>
UA2390-10 bkFusion TP Freehand	bk3000/bk5000	E14CL4b	<ul style="list-style-type: none"> <li>UA2373-03 bkFusion TP Freehand software</li> <li>UA2374 bkFusion license</li> <li>UA2371 EM sensor</li> <li>EM control unit and shelf</li> <li>EM transmitter &amp; stand</li> <li>Universal clamp</li> <li>Cables</li> </ul>

Note that sensor clamps for transducers must be ordered separately: UA2378 sensor clamp for E14CL4b

## Systems

2300-01	bk3000	2300-11	bk3000 w. battery
2300-51	bk5000	2300-61	bk5000 w. battery

2300 wheel base without battery

2300 wheel base with battery



## Upgrades

UA2608	Upgrade to SW 5.16.2 + DICOM Encrypted
UA2380	bk3000 to bk5000 Upgrade

bkSpecto = [ref] 1300

bk3000/bk5000 = [ref] 2300

bkmedical.com

BK Medical  
8 Centennial Drive  
Peabody, MA 01960  
USA  
T +1 978 326 1300  
F +1 978 326 1399

Sales and service USA  
BK Medical  
8 Centennial Drive  
Peabody, MA 01960  
USA  
T +1 978 326 1300  
F +1 978 326 1399

Europe and rest of world  
BK Medical  
Mileparken 34, 2730  
Herlev, Denmark  
T +45 4452 8100  
F +45 4452 8199



BG0565-U 2020-09

# HIGH-RESOLUTION SPECIALIZED ULTRASOUND SYSTEMS

bkSpecto | bk3000 | bk5000 | Transducers



## bkSpecto Application Packages

### UROLOGY

Application Package	Exam Types	Systems	Transducers
Small Parts Essential UA2454	Small Parts	bkSpecto	14L3, 14L3e
Prostate Essential UA2457	Prostate	bkSpecto	E11C3b, E13C2, E14CL4b
Abdomen Essential UA2458	Abdomen	bkSpecto	5C1e
Brachytherapy UA2481	Brachytherapy	bkSpecto	E14CL4b
UA2490 is required in order to enable UA2481			
Prostate Advanced UA2490	Prostate	bkSpecto	E14C4t, E10C4, E11C3b, E13C2, E14CL4b
Abdomen Advanced UA2491	Abdomen	bkSpecto	6C2, 6C2s, 9C2, 5C1e, 5P1e
Small Parts Advanced UA2492	Small Parts	bkSpecto	13L4w, 14L3, 14L3e, 18L5, X18L5s, 8L2
Essential Urology Bundle UA2494	Prostate, Small Parts, Abdomen	bkSpecto	E11C3b, E13C2, E14CL4b, 14L3, 14L3e, 5C1e
Advanced Urology Bundle UA2495	Prostate, Small Parts, Abdomen	bkSpecto	E14C4t, E10C4, E11C3b, E13C2, E14CL4b, 13L4w, 14L3, 14L3e, 18L5, X18L5s, 8L2, 6C2, 6C2s, 9C2, 5C1e, 5P1e

### PROCEDURAL ULTRASOUND

Application Package	Exam Types	Systems	Transducers
MSK & Nerve UA2453	MSK & Nerve	bkSpecto	6C2, 9C2, 5C1e, 5P1e, 13L4w, 14L3, 14L3e, 18L5, X18L5s, 8L2, 18L5s
Women's Health UA2455	OB/GYN	bkSpecto	6C2, 9C2, 5C1e, E13C2, E10C4
Vascular UA2456	Vascular	bkSpecto	9C2, 13L4w, 14L3, 14L3e, 18L5, X18L5s, 8L2, 18L5s
Colorectal UA2497	Colorectal	bkSpecto	X14L4, 20R3
Pelvic Floor UA2498	Pelvic Floor	bkSpecto	6C2, 9C2, X14L4, 20R3, E14CL4b, 5C1e, E10C4

## bkSpecto Options

UA2471	3D Freehand	UA2474	Strain Elastography	UA2479	Varian Interface
UA2472	3D Professional	UA2478	Vector Flow Imaging (VFI)		

## Systems

1300-25	bkSpecto
1300-21	bkSpecto with battery

## Upgrades

UA2496	Essential to Advanced Urology Upgrade
--------	---------------------------------------

## bk3000 / bk5000 Application Packages

### UROLOGY

Application Package	Exam Types	Systems	Transducers
Prostate UA2601	Prostate	bk3000/bk5000	E14CL4b, X14L4, E14C4t, E10C4, E13C2
Brachytherapy UA2335	Brachytherapy	bk3000/bk5000	E14CL4b
Abdomen UA2336	Abdomen	bk3000/bk5000	9C2, 6C2s, 6C2, 5P1, 5C1e
Small Parts UA2337	Small Parts	bk3000/bk5000	13L4w, 8L2, 14L3, 18L5, 18L5s, X18L5s
Urology Advanced Bundle UA2604	Prostate, Abdomen, Small Parts	bk3000/bk5000	E14CL4b, X14L4, E14C4t, E10C4, E13C2, 9C2, 6C2s, 6C2, 5P1, 5C1e, 13L4w, 8L2, 14L3, 18L5, 18L5s, X18L5s

### SURGERY

Application Package	Exam Types	Systems	Transducers
HPB Surgery UA2363	HBP Surgery	bk5000	I14C5I, I14C5T, I12C5b, I12C5f
LAP & Robotic UA2364	LAP & Robotic	bk5000	X12C4, Rob12C4k, I12C4f
Micro Neurosurgery UA2344	Micro Neuro (USA Only)	bk5000	N20P6l
Neurosurgery UA2376	Neurosurgery	bk5000	N13C5, N11C5s
Neurosurgery High Freq UA2605	Neurosurgery (USA & Canada Only)	bk5000	X18L5s <sup>d-g</sup>
Colorectal UA2339	Colorectal	bk3000/bk5000	X14L4, 20R3

UA2339 includes UA2322 3D Professional

Pelvic Floor UA2369	Pelvic Floor	bk3000/bk5000	9C2, 6C2, 5C1e, X14L4, 20R3, E14CL4b, E10C4
------------------------	--------------	---------------	---

UA2369 includes UA2322 3D Professional

### PROCEDURAL ULTRASOUND













Application Package	Exam Types	Systems	Transducers
Women's Health UA2602	Gyn/OB	bk3000/bk5000	9C2, 6C2, 6C2s, 5C1e, E14C4t, E13C2, E10C4
Vascular UA2603	Vascular	bk3000/bk5000	9C2, 5P1, 13L4w, 10L2w, 8L2, 14L3, 18L5, 18L5s, X18L5s
MSK & Nerve UA2397	MSK & Nerve	bk3000/bk5000	9C2, 6C2, 6C2s, 5P1, 5C1e, 13L4w, T7P2m <sup>h,i</sup> , 8L2, 14L3, 18L5, 18L5s, X18L5s
Pediatric UA2375	Pediatric	bk3000/bk5000	9C2, 6C2s, 8L2, 14L3, N13C5, 18L5s

Available on the bk3000 only with selection of UA2337 or UA2604

Available on the bk5000 only with selection of UA2376








## bk3000 / bk5000 Options

UA2321	3D Freehand	UA2327	Varian Interface	UA2328	Elastography
UA2322	3D Professional	UA2333	Needle Enhancement (X-Shine)	UA2331	WiFi
UA2324	Contrast Enhanced	UA2334	Brainlab Neuro Navigation	UA2374	bkFusion License
UA2325	Vector Flow Imaging (VFI)	UA2340	Medtronic Liver Navigation	UA2389	DICOM Encrypted

	Prostate Triplane E14C4t (9018)	High-Resolution Endfire E13C2 (9029)	Endocavity Biplane E14CL4b (9048)	High-Resolution Biplane E11C3b (9008)	Endocavity E10C4 (9019)	Endocavity 3D X14L4 (9038)	Anorectal 3D 20R3 (9052)	Linear 14L3 (9051)	Linear 8L2 (9032)	Wide Linear 10L2w (9022)	High-Frequency Linear 18L5 (9070)
											
MAIN USES (System Dependent)	Transrectal, Transvaginal	Transrectal, Transvaginal	Transrectal, Transvaginal	Transrectal	Transvaginal, Transrectal, Fetal incl. OB	Transrectal, Transvaginal	Transrectal, Transvaginal	Small Parts, MSK, PV, Pediatrics	PV, Pediatrics, Small Parts, MSK	PV	Small Parts, MSK, PV
Puncture	•	•	•	•	•			•	•	•	•
Contrast Imaging	•	•	•			•			•		
Smart Button™	•	•	•	•		•	•	•	•	•	•
Elastography	•	•	•	•	•	•		•	•	•	•
VFI								•	•	•	
Tru-Focus									•		
Sterilizable*	•	•	•	•	•	• <sup>a</sup>	• <sup>a</sup>	•	•	•	•











OPERATES WITH FOLLOWING SYSTEMS:

bk5000	•	•	•		•	•	•	•	•	•	•
bk3000	•	•	•		•	•	•	•	•	•	•
bkSpecto	•	•	•	•	•	•	•	•	•		•

	Small High-Frequency Linear 18L5s (9081)	Linear 14L3e (9086)	Wide Linear 13L4w (9011)	High-Frequency Curved 9C2 (9002)	Curved 5C1e (9085)	Curved 6C2 (9040)	Small Curved 6C2s (9023)	Small Footprint Cardiac 5P1 (9077)	Small Footprint Cardiac 5P1e (9087)	Hockey Stick X18L5s (9009)	TEE T7P2m (9027) <sup>h,i</sup>
											
MAIN USES (System Dependent)	MSK, Small Parts, PV, Pediatrics	Small Parts, MSK, Pediatrics, PV	Small Parts, MSK, PV, Pediatrics	Abdominal, MSK, Fetal incl. OB, PV, Pediatrics	Abdominal, MSK, Fetal incl. OB	Abdominal, MSK, Fetal incl. OB	Abdominal, Pediatrics, Fetal incl. OB	Cardiac <sup>b</sup> , Abdominal	Cardiac <sup>b</sup> , Abdominal, Pediatric	Intraoperative <sup>g</sup> , Intraoperative (neuro) <sup>g</sup> , PV, MSK, Small Parts	Transesophageal Cardiac
Puncture		•	•	•	•	•	•				
Contrast Imaging						•	•				
Smart Button™	•	•	•	•	•	•	•	•	•	•	
Elastography	•	•	•	•						•	
VFI		•									
Tru-Focus					•	•					
Sterilizable*			•	•	•	•	•	•	•	•	

OPERATES WITH FOLLOWING SYSTEMS:

bk5000	•		•	•	•	•	•	•	•	•	• <sup>e,h,i</sup>
bk3000	•		•	•	•	•	•	•	•	•	• <sup>e,j</sup>
bkSpecto	•	•	•	•	•	•	•		•	•	

	4-Way Laparoscopic I12C4f (9066)	Intraoperative I14C5I (9015)	T-shaped Intraoperative I14C5T (9016)	Intraoperative Biplane I12C5b (9024)	Mini-T Intraoperative I12C5 (9034)	Robotic Drop-In X12C4 (9026)	Robotic Drop-In Rob12C4 (9096) <sup>k</sup>	Craniotomy N13C5 (9062)	Burr Hole N11C5s (9063)	Minimally Invasive N20P6 (9007)
										
MAIN USES (System Dependent)	Intraoperative	Intraoperative, Pediatrics	Intraoperative, Abdominal	Intraoperative, MSK, PV	Intraoperative, MSK, Small Parts	Intraoperative	Intraoperative	Intraoperative, Intraoperative (neuro), Pediatric Neonatal Cephalic	Intraoperative, Intraoperative (neuro),	Intraoperative (neuro)
Puncture	•	•		•				•	•	
Contrast Imaging	•	•	•	•	•	•	•			
Smart Button™								•	•	•
Elastography <sup>c</sup>	•	•	•	•	•	•	•	•		
VFI										
Tru-Focus	•	•	•			•	•			
Sterilizable*	•	•	•	•	•	•	•	•	•	•

OPERATES WITH FOLLOWING SYSTEMS:

bk5000	•	•	•	•	• <sup>e</sup>	•	• <sup>e</sup>	•	•	• <sup>e</sup>
bk3000								•		
bkSpecto										

MSK = Musculoskeletal PV = Peripheral Vascular OB = Obstetrics  
 \* For complete information about sterilization compatibility and validity, please see Care & Cleaning.  
 b: Not for direct use on the heart.  
 a: Sterilization with STERIS System 1 models only.  
 c: Use of N13C5 and X18L5s for elastography has not been CE-marked.  
 d: Use of X18L5s for intraoperative (neuro) has not been CE-marked.  
 e: This transducer is compatible with software version 5.16.2 and higher.  
 f: Only available in USA.

g: Use of X18L5s for intraoperative has not been cleared by Health Canada for use on common carotid  
 h: T7P2m has not been licensed by Health Canada for use on the bk5000  
 i: T7P2m is not CE-marked.  
 j: I12C5 has not been licensed by Health Canada.  
 k: Rob12C4 has not been licensed by Health Canada.  
 l: N20P6 has not been licensed by Health Canada and is not CE-marked.