



Premium Performance – Exceptional Value





Powerful, Portable, Complete – the Best Fit for Your Busy Clinic

- High Resolution Imaging
- Advanced Imaging Technologies
- · Streamlined Workflow and Improved Productivity
- Comprehensive Range of Innovative Transducers

Maximize Your Imaging Performance

Full-featured ultrasound system designed for continuous imaging performance, mobility and ease of use

High Resolution Imaging

- Premium performance with efficiency and speed
- High resolution 19" monitor
- Sensitive Color Doppler with superb spatial resolution and sensitivity

Advanced Imaging Technologies

- Motion Compensated Angular Compound Imaging (MACI) ensures high quality real-time images using camera imaging technology
- Automatic Mode Adjustment (AMA) optimizes high resolution B-Mode imaging when you change modes and depth
- Triplex enables you to maximize PRF with a tight spectral envelope and real-time synchronicity between 2D and Color Doppler

Intuitive and Mobile

- Small footprint fits in the tightest spaces
- Four hours plug-free imaging
- Intuitive and easy to use
- Smooth, sealed keyboard for easy cleaning and disinfection

State-of-the-Art Transducers for Your Clinical Needs

- Innovative design for access to all areas
- Advanced puncture guides
- Convenient one-button control
- Easy sterilization and disinfection



Urology

Gold standard for all urological imaging. Best solution for prostate, kidney and testicular imaging.

The unique triplane transducer with simultaneous biplane and endfire imaging capabilities is ideal for quiding biopsies.



Colorectal

See every detail.

View all layers of the rectal wall.

Get detailed views of anal sphincter and other clinically relevant regions with 2D and 3D images.




Anesthesia/Musculoskeletal

Superior nerve and needle visualization.

Ideal for regional anesthesiology, pain clinics and MSK applications.

Dedicated application-specific interface.



Vascular

Unique Vector Flow Imaging mode enables angle-independent visualization of blood flow.



Pelvic Floor

Get a complete picture of the entire pelvic floor in 2D and 3D.

High resolution vaginal, anal and transperineal ultrasound imaging – quick and easy.




General Imaging

Exceptional image quality, intuitive user interface and advanced workflow features enable quality abdominal, small part and vascular imaging.

Real Clinical Impact with High Resolution Imaging




Hydrocele and Testis




Transverse Prostate with Color Doppler Imaging



Kidney with Solid Mass on Lower Pole



PW Doppler of SMA using Triplex Mode




PW Doppler of Distal CCA



Solid Renal Mass on Lower Pole




PW Doppler of Right Renal Artery



Biometry Measurements of Prostate using Simultaneous Biplane Imaging




Color Doppler Imaging of Abdominal Aorta, Celiac Trunk and SMA




Rectal Wall Tumor




Color Doppler Imaging of Hepatic Vessels



Biopsy-Proven Malignant Breast Mass




Baker's Cyst



3D Volume Acquisition of Prostate



Sagittal Liver and IVC



PW Doppler of Renal Arcuate Artery

Acoustic and Clinical Innovators for Over 40 Years

Recognized as a global leader in the development of ultrasound systems designed for surgery, urology and anesthesia, BK Ultrasound offers a versatile range of state-of-the art imaging systems, technologies and dedicated transducers to help enhance diagnostic confidence and advance the care of patients.

Prostate Triplane	Prostate Biplane	3D ART	Endocavity Biplane	Endfire Curved	Curved
8818	8808e	8838	8848	8667	8823
12–4 MHz	10-5 MHz	12–4 MHz	12–4 MHz	10–5 MHz	6–1.8 MHz
Curved	Curved	Linear	Linear	High Frequency	Vascular
8820e	8830	8811	8670	Linear 8870	8822
6–2 MHz	6–2 MHz	12–5 MHz	12–4 MHz	18–6 MHz	9–3.5 MHz
	-		-		1
Anorectal 3D	Rigid Laparoscopic	Craniotomy	"Hockey Stick"	Small footprint	Endovaginal
2052	8836*	8862	8809	Cardiac 8837	8819
2052 16–6 MHz	8836* 12-5 MHz	8862 10–3.8 MHz	8809 15–6 MHz		8819 9–5 MHz

^{*} In selected countries only.

System versions available: flex Focus 400, flex Focus 400 exp, flex Focus 500

Direct Sales Offices

BK Ultrasound				
8 Centennial Drive				
Peabody				
MA 01960				
USA				
T+1 978-326-1300				
F+1978-326-1399				
bkultrasound.com				

Ana	logic	Bene	lux	N.V

Generaal de Wittelaan 19B bus 7 2800 Mechelen Belaium T+32 15 280 600 F +32 15 291 025 bkultrasound.com

Mileparken 34 2730 Herlev Denmark T +45 4452 8100 F +45 4452 8199 bkultrasound.com

BK Medical ApS

25451 Quickborn Germany bkultrasound.com

bk medical Medizinische Systeme GmbH

Pascalkehre 13 T +49 4106 99 55-0 F +49 4106 99 55-99

Analogic Italia S.r.l.

Via Morandi, 10 20077 Melegnano MI T +39 02 90 78 13 47 F +39 02 90 78 19 05 bkultrasound.com

B-K Medical AB

Västberga Allé 26, 3tr 126 30 Hägersten Sweden T +46 8-744 02 11 F +46 8-744 02 12 bkultrasound.com

BK Medical

Unit A20/A20a Basepoint Business & Innovation Centre 110 Butterfield, Great Marlings Luton Bedfordshire LU2 8DL United Kingdom T+44 1582 433710 F +44 1582 433716 bkultrasound.com



BK Medical ApS, Mileparken 34, 2730 Herlev, Denmark. T +45 4452 8100 F +45 4452 8199

BK Ultrasound 8 Centennial Drive Peabody MA 01960 T+1 978-326-1300 bkultrasound.com **North America** Sales & Service **BK Ultrasound** 8 Centennial Drive Peabody

USA T+1 978-326-1300 F+1978-326-1399 bkultrasound.com

Europe and Rest of World Sales, Service & Design Center

BK Ultrasound Mileparken 34 2730 Herlev Denmark T+45 4452 8100 F +45 4452 8199 bkultrasound.com

Asia Sales & Service

Analogic Medical Equipment (Shanghai) Co., Ltd. 1377, Lan Dian Road Pu Dong New District Shanghai China 201132 T+86 21 2089 0333 bkultrasound.com

powered by analogic 🚇

Canada **Design Center**

Ultrasonix Medical Corporation 130-4311 Viking Way Richmond BC V6V 2K9 Canada T+1866-437-9508 F+1604-279-8559 bkultrasound.com

MA 01960