

A New Level of Ultrasound Guidance for Your Procedures



- Upper Extremity Nerve Blocks
- Lower Extremity Nerve Blocks
- Spinal Nerve Blocks
- TAP Blocks
- Vascular Access




bkSpecto

Anesthesiology & Pain Medicine


40 Years of Expertise in Specialized Ultrasound Guidance

Outstanding Image Quality to Guide Your Procedures


Upper Extremity Blocks




Ulnaris nerve - needle
visualization
18L5 transducer



Median nerve with color Doppler
18L5 transducer



Interscalene
18L5 transducer



Brachial plexus supraclavicular
13L4w transducer



14L3/e (9051/9086)
Linear



13L4w (9011)
Wide Linear




18L5 (9070)
High-Frequency
Linear




8L2 (9032)
Linear


Lower Extremity Blocks




Popliteal nerve
with color Doppler
18L5 transducer



Sciatic nerve
9C2 transducer



Popliteal nerve
18L5 transducer



Saphenous nerve
9C2 transducer



14L3/e (9051/9086)
Linear



13L4w (9011)
Wide Linear



18L5 (9070)
High-Frequency
Linear



8L2 (9032)
Linear




9C2 (9002)
Curved




5C1e (9085)
Curved

Spinal/TAP Blocks



Spine
9C2 transducer



Muscular layers of the anterior
abdominal wall (TAP)
18L5 transducer



9C2 (9002)
Curved



5C1e (9085)
Curved



8L2 (9032)
Linear




18L5 (9070)
High-Frequency
Linear




14L3/e (9051/9086)
Linear

Vascular Access



Radial artery with color
Doppler
18L5s transducer




Proximal ICA PW Doppler
8L2 transducer



14L3/e (9051/9086)
Linear



8L2 (9032)
Linear



18L5s (9081)
Linear



X18L5s (9009)
Hockey Stick




18L5 (9070)
High-Frequency
Linear

Specialty Transducers



Median nerve
18L5s transducer



Radial artery and nerve
X18L5s transducer



X18L5s (9009)
Hockey Stick



18L5s (9081)
Linear



5P1e (9087)*
Small Footprint
Cardiac

Fully flexible
transducer tip for
convenient
imaging,
including for
vascular access.

High-frequency,
small footprint
transducer for
superficial nerve
blocks and
pediatrics.

Ideal transducer
for cardiac
imaging and
FAST/FATE
exams.

* 5P1e is not for direct use on the heart

Personalize Your Workflow and Simplify Procedures

bkSpecto


- Excellent image quality and needle visualization for fast and efficient procedures.
- Access to a range of BK's high-performance, sterilizable transducers for advanced procedures.
- Integrated battery and ergonomic design for mobility in your procedure room.
- Access to advanced technology and performance modes: Vector Flow Imaging (VFI), Needle Enhance, DICOM, and WIFI.



Personalize the unique interface to perform your blocks and vascular access procedures with ease.



Sealed glass interface makes cleaning and infection control fast and efficient.



Custom-fit clinical design. Adjust monitor orientation to suit your needs.

BK Medical
8 Centennial Drive
Peabody
MA 01960
USA
T +1 978-326-1300
bkmedical.com

**USA
Sales & Service**
BK Medical
8 Centennial Drive
Peabody
MA 01960
USA
T +1 978-326-1300
F +1 978-326-1399
bkmedical.com

**Europe and Rest of World
Sales, Service & Design Center**
BK Medical
Mileparken 34
2730 Herlev
Denmark
T +45 4452 8100
F +45 4452 8199
bkmedical.com