

ULTRASOUND DESIGNED

bk5000


FOR SURGERY

- Excellent images instantly - no adjustment needed
- Small system footprint - fits easily in the OR
- Wide range of specialized and sterilizable transducers
- Sterilizable remote control for convenient control of the system in the sterile field
- Rapid boot time and up to two hours of plug-free imaging
- Advanced modes: elastography, contrast imaging, 3D, vector flow imaging (VFI)



bk 
medical


HPB, INTRAOPERATIVE



Thrombus in portal vein
I12C5b transducer

Elastography of metastatic lesion in liver
I14C5T transducer

I14C5i (9015) Intraoperative
I14C5T (9016) T-Shaped Intraoperative
I12C5b (9024) Intraoperative Biplane
I12C5 (9034)¹ Mini-T Intraoperative
I12C4f (9066) 4-Way Laparoscopic



NEUROSURGERY, INTRAOPERATIVE



Bilateral ventriculomegaly
N11C5s transducer

Metastatic lesion
N13C5 transducer

N13C5 (9062) Craniotomy
N11C5s (9063) Burr Hole
X18L5s (9009)² Hockey Stick
N20P6 (9007)^{3,4} Minimally Invasive



ROBOTIC-ASSISTED




Renal tumor
Rob12C4 transducer

Rob12C4 (9096)¹ Robotic




ENDOCAVITY




Expanded sector imaging of prostate gland
E14C4t transducer

bkFusion Biopsy of region of interest
E14C4t transducer

E14C4t (9018) Prostate Triplane
E13C2 (9029) High Resolution Endfire
E14CL4b (9048) Endocavity Biplane
E10C4 (9019) Endocavity




ANORECTAL, PELVIC FLOOR



Urethral mass
X14L4 transducer

9C2 (9002) High-Frequency Curved
E14CL4b (9048) Endocavity Biplane
X14L4 (9038) Endocavity 3D
20R3 (9052) Anorectal 3D



LAPAROSCOPIC



Color Doppler of liver hemangioma
I12C4f transducer

I12C4f (9066) 4-Way Laparoscopic



PHASED

5P1 (9077)⁴ Small Footprint Cardiac
T7P2m (9027)¹ TEE



CURVED

6C2 (9040) Curved
5C1e (9085) Curved
6C2s (9023) Small Curved
9C2 (9002) High Frequency Curved



LINEAR

13L4w (9011) Wide Linear
14L3 (9051) Linear
10L2w (9022) Wide Linear
8L2 (9032) Linear
18L5 (9070) High Frequency Linear
18L5s (9081) Small High Frequency Linear



¹This transducer is compatible with bk5000 software version 5.16.2 and higher; This transducer has not been licensed by Health Canada for use on the bk5000.

²The Hockey Stick Transducer X18L5s is not approved for use on carotid in Canada; Use of the Hockey Stick transducer for intraoperative (neuro) has not been CE-marked.

³The Minimally Invasive Transducer N20P6 has not been licensed by Health Canada or CE-marked.

⁴The Cardiac Transducer 5P1 is not for direct use on the heart.

BK Medical
8 Centennial Drive
Peabody,
MA 01960
USA
T +1 978 326 1300
bkmedical.com

USA
Sales & Service
BK Medical
8 Centennial Drive
Peabody, MA 01960
USA
T +1 978 326 1300
F +1 978 326 1399
bkmedical.com

Europe and Rest of World
Sales, Service & Design Center
BK Medical
Mileparken 34,
2730 Herlev, Denmark
T +45 4452 8100
F +45 4452 8199
bkmedical.com

bk
medical