

# CHOOSE THE PACKAGE THAT FITS YOUR NEEDS

SOLUTION	ORDER NUMBER	INCLUDED
bk5000 Hepatobiliary with battery	bk5000b-HPB	bk5000 (2300-61), I12C4f, I14C5T, I12C5b, UA2323, UA2302
bk5000 Hepatobiliary	bk5000-HPB	bk5000 (2300-51), I12C4f, I14C5T, I12C5b, UA2323, UA2302
bk5000 Neurosurgery with battery	bk5000b-NEURO	bk5000 (2300-61), N13C5, N11C5s, UA2323, UA2302
bk5000 Neurosurgery	bk5000-NEURO	bk5000 (2300-51), N13C5, N11C5s, UA2323, UA2302
bk5000 Robotic-Assisted Surgery with battery	bk5000b-RAS	bk5000 (2300-61), X12C4, UA2323, UA2302
bk5000 Robotic-Assisted Surgery	bk5000-RAS	bk5000 (2300-51), X12C4, UA2323, UA2302
bk3000 Urology	bk3000-URO	bk3000 (2300-01), E14C4t, 6C2, 13L4w, UA2301, UA2323, UA2360
bk3000 Urology Interventional	bk3000-URO-INT	bk3000 (2300-01), E14C4t, 6C2s, 13L4w, UA2301, UA2323, UA2360

## THANK YOU

for considering **bk** – where we excel at tailoring the best ultrasound imaging for your unique needs.

[bkultrasound.com](http://bkultrasound.com)



BK Medical  
8 Centennial Drive  
Peabody, MA 01960  
USA  
T +1 978 326 1300

Sales and service USA  
BK Medical  
8 Centennial Drive  
Peabody, MA 01960  
USA  
T +1 978 326 1300  
F +1 978 326 1399

Europe and rest of world  
BK Medical  
Mileparken 34, 2730  
Herlev, Denmark  
T +45 4452 8100  
F +45 4452 8199

BK Medical Aps, Mileparken 34, 2730 Herlev, Denmark. T +45 4452 8100 F +45 4452 8199

BG0565-J 09/18

# HIGH RESOLUTION SPECIALIZED ULTRASOUND SYSTEMS

bkSpecto  
bk3000  
bk5000  
Transducers



# APPLICATION PACKAGES

SURGERY			
Application package	Pro Packages/ Exam Types	Systems	Transducers
Surgery UA2302	Neuro	bk5000	N13C5, N11C5s
	Surgery	bk5000	14C5I, I14C5T, I12C5b, I12C4f, X12C4, X18L5s
UA2302: Enabling Intraoperative scanning. HPB, Neuro and Robotic Assisted Surgery			
Partial Nephrectomy UA2312	Surgery	bk5000	X12C4, I12C4f
UA2312: Subset of UA2302 enabling Robotic Assisted Surgery			
Colorectal & Pelvic Floor UA2308	Colorectal	bk3000/5000	20R3, X14L4
	Pelvic Floor	bk3000/5000	20R3, X14L4, E14CL4b, E10C4, E13C2, 9C2, 6C2, 6C2s, 5C1e
UA2308: Enabling Colorectal and Pelvic Floor 3D diagnostics. Includes UA2322 3D Professional			
SUR Super AppPack UA2316	Includes the following application packages	UA2301 Urology, UA2302 Surgery, UA2303 Breast/Prostate Therapy, UA2304 General Imaging, UA2305 Women's Health, UA2306 Peripheral Vascular, UA2309 Point of Care, UA2311 Emergency Medicine	
UROLOGY			
Application package	Pro Packages/ Exam Types	Systems	Transducers
Urology UA2301	URO Prostate	bk3000/5000	E14C4t, E14C4, E13C2, E14CL4b, E10C4
	URO Abdominal	bk3000/5000	6C2, 6C2s, 5C1e
	URO Small Parts	bk3000/5000	13L4w, 14L3, 18L5, 18L5s, 8L2
Breast/Prostate Therapy UA2303	Brachytherapy	bk3000/5000	E14CL4b
	URO Prostate	bk3000/5000	E14C4t, E13C2
	Breast	bk3000/5000	13L4w, 18L5, 14L3, 8L2
URO Super AppPack UA2315	Includes the following application packages	UA2301 Urology, UA2303 Breast/Prostate Therapy, UA2304 General Imaging, UA2305 Women's Health, UA2306 Peripheral Vascular, UA2309 Point of Care, UA2311 Emergency Medicine	
Urology Basic UA2450	Prostate, Abdomen, Small Parts	bkSpecto	5C1e, 14L3e, E11C3b, E13C2
Urology Premium UA2451	Prostate, Abdomen, Small Parts	bkSpecto	5C1e, 18L5, 14L3e, E14C4t, E11C3b, E13C2, E14CL4b, E10C4, 6C2, 6C2s, 9C2, 14L3, 13L4w

POINT OF CARE			
Application package	Pro Packages/ Exam Types	Systems	Transducers
Point of Care UA2309	UGRA	bk3000/5000	6C2s, 9C2, 13L4w, 18L5, 14L3, 10L2w, 6C2, 5C1e, 8L2, 18L5s
	MSK	bk3000/5000	6C2s, 9C2, 13L4w, 18L5, 14L3, 10L2w, 6C2, 8L2, 5C1e, 18L5s, X18L5s
	Cardiac	bk3000/5000	5P1
	GYN	bk3000/5000	E10C4, 6C2, 6C2s, 9C2, E14C4t, E14C4, E13C2, 5C1e
Emergency Medicine UA2311	EMED Abdomen	bk3000/5000	6C2, 6C2s, 5C1e, 9C2, 5P1
	EMED Vascular	bk3000/5000	18L5, 13L4w, 14L3, 10L2w, 8L2, 9C2, 5P1
	EMED Small Parts	bk3000/5000	13L4w, 18L5, 14L3, 8L2, 18L5s
	EMED Cardiac	bk3000/5000	5P1
	EMED OB	bk3000/5000	9C2, E10C4, 6C2, 5C1e
SHARED SYSTEMS			
Application package	Pro Packages/ Exam Types	Systems	Transducers
General Imaging UA2304	Abdominal	bk3000/5000	6C2, 6C2s, 5C1e, 9C2, 5P1
	Small Parts	bk3000/5000	13L4w, 18L5, 14L3, 8L2, 18L5s
	Pediatrics	bk3000/5000	6C2s, 9C2, 5C1e, 18L5, 18L5s, N13C5
Women's Health UA2305	Breast	bk3000/5000	13L4w, 18L5, 14L3, 8L2
	OB	bk3000/5000	E10C4, 6C2, 9C2, 5C1e
	GYN	bk3000/5000	E10C4, 6C2, 6C2s, 9C2, E14C4t, E14C4, E13C2, 5C1e
	GYN Follicle	bk3000/5000	E10C4, E13C2, 6C2, 6C2s, 9C2, 5C1e
UA2306 Peripheral Vascular	Peripheral Vascular	bk3000/5000	9C2, 14L3, 13L4w, 10L2w, 18L5, 8L2, 18L5s
	Carotid	bk3000/5000	18L5, 13L4w, 14L3, 10L2w, 8L2

	Prostate Triplane E14C4t (9018)	High Resolution Endfire E13C2 (9029)	Endfire Curved E14C4 (9067)	Endocavity Biplane E14CL4b (9048)	High Resolution Biplane E11C3b (9008)	Endocavity E10C4 (9019)	Endocavity 3D X14L4 (9038)	Anorectal 3D 20R3 (9052)	Linear 14L3 (9051)	Linear 8L2 (9032)
										
MAIN USES (System Dependent)	Transrectal, Transvaginal	Transrectal, Transvaginal	Transrectal, Transvaginal	Transrectal, Transvaginal	Transrectal	Transvaginal, Transrectal, Fetal incl. Obstetrics	Transrectal, Transvaginal	Transrectal, Transvaginal	Small Parts, MSK, PV, Pediatrics	PV, Pediatrics, Small Parts, MSK
Footprint/Tip diameter (mm)	20	20	22	20	20	22 x 18	20	17	45 x 14	57 x 10
Puncture	•	•	•	•	•	•	•		•	•
Contrast Imaging <sup>a)</sup>	•			•			•			•
Smart Button™	•	•	•	•	•		•	•	•	•
Elastography	•	•	•	•	•		•		•	•
VFI <sup>b)</sup>									•	•
Sterilizable**	•	•	•	•	•	•	• <sup>e)</sup>	• <sup>e)</sup>	•	•


OPERATES WITH FOLLOWING SYSTEMS

bkSpecto*	•	•		•	•	•			•	
bk3000	•	•	•	•		•	•	•	•	•
bk5000	•	•	•	•		•	•	•	•	•

	High Frequency Linear 18L5 (9070)	Small High Frequency Linear 18L5s (9081)	Wide Linear 10L2w (9022)	Linear 14L3e (9086)	Wide Linear 13L4w (9011)	High Frequency Curved 9C2 (9002)	Curved 5C1e (9085)	Curved 6C2 (9040)	Small Curved 6C2s (9023)
									
MAIN USES (System Dependent)	Small Parts, MSK, PV	MSK, Small Parts, PV, Pediatrics	PV	Small Parts, MSK, Pediatrics, PV	Small Parts, MSK, PV, Pediatrics	Abdominal, MSK, Fetal incl. Obstetrics, PV, Pediatrics	Abdominal, Obstetrics, Fetal, MSK	Abdominal, MSK, Fetal incl. Obstetrics	Abdominal, Pediatrics, Fetal incl. Obstetrics
Footprint/Tip diameter (mm)	48 x 13	24 x 11	57 x 10	45 x 14	57 x 10	54 x 10	68 x 19	69 x 19	33 x 12
Puncture	•		•	•	•	•		•	•
Contrast Imaging <sup>a)</sup>							•	•	•
Smart Button™	•	•	•	•	•	•	•	•	•
Elastography <sup>d)</sup>	•	•	•	•	•	•			
VFI <sup>b)</sup>			•						
Sterilizable**	•		•		•	•		•	•

OPERATES WITH FOLLOWING SYSTEMS

bkSpecto*	•			•	•	•	•	•	•
bk3000	•	•	•		•	•	•	•	•
bk5000	•	•	•		•	•	•	•	•

	4-Way Laparoscopic I12C4f (9066)	Intraoperative I14C5I (9015)	T-Shaped Intraoperative I14C5T (9016)	Intraoperative Biplane I12C5b (9024)	Robotic Drop-In X12C4 (9026)	Craniotomy N13C5 (9062)	Burr-Hole N11C5s (9063)	"Hockey-Stick" X18L5s (9009)	Small Footprint Cardiac 5P1 (9077)
									
MAIN USES (System Dependent)	Intraoperative	Intraoperative, Pediatrics	Intraoperative, Abdominal	Intraoperative, MSK, PV	Intraoperative	Intraoperative, Neurosurgery, Neonatal Cephalic, Pediatrics	Intraoperative, Neurosurgery	Intraoperative, PV, Small Parts, MSK	Cardiac <sup>d)</sup> , Abdominal, Transcranial
Footprint/Tip diameter (mm)	33 x 9	60 x 14	60 x 14	40 x 30	33 x 9	29 x 10	10 x 9	31 x 9	26 x 20
Puncture	•	•	•	•	•	•	•	•	
Contrast Imaging <sup>a)</sup>	•	•	•	•	•				
Smart Button™						•	•	•	•
Elastography	•	•	•	•	•	•	•	•	
VFI <sup>b)</sup>									
Sterilizable**	•	•	•	•	•	•	•	•	•

OPERATES WITH FOLLOWING SYSTEMS

bkSpecto*									
bk3000						•		•	•
bk5000	•	•	•	•	•	•	•	•	•

a) Contrast imaging on bk3000/bk5000 has not been licensed by Health Canada.  
b) VFI on bk3000/bk5000 has not been licensed by Health Canada.  
c) Not for direct use on the heart.  
d) Use of the N13C5, N11C5s and X18L5s for elastography has not been CE-marked.  
e) Sterilization with STERIS System 1 models only.

MSK = Musculoskeletal PV = Peripheral Vascular  
Health Canada: E11C3b, 14L3e not licensed. I12C4f, I14C5I, I14C5T, X12C4, I12C5b only licensed for BK5000. E14C4, 9C2, 8L2, 20R3, X14L4, X18L5s, 5C1e, E13C2 licensed for bk3000/5000.

\* bkSpecto is not licensed by Health Canada  
\*\* For complete information about sterilization compatibility and validity, please see Care, Cleaning & Safety.