

APPLICATIONS PACKAGES CONTINUED

SHARED SYSTEMS

Application package	Pro Packages/ Exam Types	Systems	Transducers
General Imaging UA2304	Abdominal	bk3000/5000	6C2, 6C2s, 5C1e, 9C2, 5P1
	Small Parts	bk3000/5000	13L4w, 18L5, 14L3, 8L2, 18L5s
	Pediatrics	bk3000/5000	6C2s, 9C2, 18L5, 18L5s, N13C5
Women's Health UA2305	Breast	bk3000/5000	13L4w, 18L5, 14L3, 8L2
	OB	bk3000/5000	E10C4, 6C2, 9C2, 5C1e
	GYN	bk3000/5000	E10C4, 6C2, 6C2s, 9C2, E14C4t, E13C2, 5C1e
	GYN Follicle	bk3000/5000	E10C4, E13C2, 6C2, 6C2s, 9C2, 5C1e
Peripheral Vascular UA2306	Peripheral Vascular	bk3000/5000	9C2, 14L3, 13L4w, 18L5, 8L2, 18L5s, 10L2w
	Carotid	bk3000/5000	18L5, 13L4w, 14L3, 8L2, 10L2w

CHOOSE THE PACKAGE THAT FITS YOUR NEEDS

SOLUTION	ORDER NUMBER	INCLUDED
bk5000 Hepatobiliary with battery	bk5000b-HPB	bk5000 (2300-61), I12C4f, I14C5T, I12C5b, UA2323, UA2302
bk5000 Hepatobiliary	bk5000-HPB	bk5000 (2300-51), I12C4f, I14C5T, I12C5b, UA2323, UA2302
bk5000 Neurosurgery with battery	bk5000b-NEURO	bk5000 (2300-61), N13C5, N11C5s, UA2323, UA2302
bk5000 Neurosurgery	bk5000-NEURO	bk5000 (2300-51), N13C5, N11C5s, UA2323, UA2302
bk5000 Robotic-Assisted Surgery with battery	bk5000b-RAS	bk5000 (2300-61), X12C4, UA2323, UA2302
bk5000 Robotic-Assisted Surgery	bk5000-RAS	bk5000 (2300-51), X12C4, UA2323, UA2302
bk3000 Urology Urology	bk3000-URO	bk3000 (2300-01), E14C4t, 6C2, 13L4w, UA2301, UA2323, UA2360
bk3000 Urology Interventional	bk3000-URO-INT	bk3000 (2300-01), E14C4t, 6C2s, 13L4w, UA2301, UA2323, UA2360

THANK YOU
for considering bk – where we excel at tailoring the best ultrasound imaging for your unique needs.

bkmedical.com

BK Medical
8 Centennial Drive
Peabody, MA 01960
USA
T +1 978 326 1300
F +1 978 326 1399

Sales and service USA
BK Medical
8 Centennial Drive
Peabody, MA 01960
USA
T +1 978 326 1300
F +1 978 326 1399

Europe and rest of world
BK Medical
Mileparken 34, 2730
Herlev, Denmark
T +45 4452 8100
F +45 4452 8199



BG0565-M 05/19

HIGH-RESOLUTION SPECIALIZED ULTRASOUND SYSTEMS

bk5000
bk3000
bkSpecto
Transducers



APPLICATION PACKAGES

SURGERY

Application package	Pro Packages/ Exam Types	Systems	Transducers
UA2302: Enabling Intraoperative scanning. HPB, Neuro and Robotic Assisted Surgery			
Surgery UA2302	Neuro	bk5000	N13C5, N11C5s
	Surgery	bk5000	I14C5i, I14C5T, I12C5b, I12C4f, X12C4, X18L5s
UA2312: Subset of UA2302 enabling Robotic Assisted Surgery			
Partial Nephrectomy UA2312	Surgery	bk5000	X12C4, I12C4f
UA2308: Enabling Colorectal and Pelvic Floor 3D diagnostics. Includes UA2322 3D Professional			
Colorectal & Pelvic Floor UA2308	Colorectal	bk3000/5000	20R3, X14L4
	Pelvic Floor	bk3000/5000	20R3, X14L4, E14CL4b, E10C4, E13C2, 9C2, 6C2, 6C2s, 5C1e
SUR Super AppPack UA2316	Includes the following application packages	UA2301 Urology, UA2302 Surgery, UA2303 Breast/Prostate Therapy, UA2304 General Imaging, UA2305 Women's Health, UA2306 Peripheral Vascular, UA2309 Point of Care, UA2311 Emergency Medicine	
Brain and Spinal Cord UA2317	Spinal Cord	bk5000	X18L5s ^o

UROLOGY











Application package	Pro Packages/ Exam Types	Systems	Transducers
Urology UA2301	URO Prostate	bk3000/5000	E14C4t, E13C2, E14CL4b, E10C4
	URO Abdominal	bk3000/5000	6C2, 6C2s, 5C1e
	URO Small Parts	bk3000/5000	13L4w, 14L3, 18L5, 8L2
Breast/Prostate Therapy UA2303	Brachytherapy	bk3000/5000	E14CL4b
	URO Prostate	bk3000/5000	E14C4t, E13C2
	Breast	bk3000/5000	13L4w, 18L5, 14L3, 8L2
URO Super AppPack UA2315	Includes the following application packages	UA2301 Urology, UA2303 Breast/Prostate Therapy, UA2304 General Imaging, UA2305 Women's Health, UA2306 Peripheral Vascular, UA2309 Point of Care, UA2311 Emergency Medicine	
Small Parts Essential UA2454	Small Parts	bkSpecto*	14L3e
Small Parts Advanced UA2492	Small Parts	bkSpecto*	13L4w, 14L3, 14L3e, 18L5, X18L5s, 8L2

UROLOGY CONTINUED

Application package	Pro Packages/ Exam Types	Systems	Transducers
Prostate Essential UA2457	Prostate	bkSpecto*	E11C3b, E13C2, E14CL4b
Prostate Advanced UA2490	Prostate	bkSpecto*	E14C4t, E10C4, E11C3b, E13C2, E14CL4b
Abdomen Essential UA2458	Abdomen	bkSpecto*	5C1e
Abdomen Advanced UA2491	Abdomen	bkSpecto*	6C2, 6C2s, 9C2, 5C1e, 5P1e
Brachytherapy UA2481	Brachytherapy	bkSpecto*	E14CL4b

PROCEDURAL ULTRASOUND

Application package	Pro Packages/ Exam Types	Systems	Transducers
Point of Care UA2309	UGRA	bk3000/5000	6C2s, 9C2, 13L4w, 18L5, 14L3, 6C2, 5C1e, 8L2, 18L5s, 10L2w
	MSK	bk3000/5000	9C2, 13L4w, 18L5, 14L3, 6C2, 8L2, 5C1e, 18L5s, X18L5s, 10L2w
	Cardiac	bk3000/5000	5P1
	GYN	bk3000/5000	E10C4, 6C2, 6C2s, 9C2, E14C4t, E13C2, 5C1e
Emergency Medicine UA2311	EMED Abdomen	bk3000/5000	6C2, 6C2s, 5C1e, 9C2, 5P1
	EMED Vascular	bk3000/5000	18L5, 13L4w, 14L3, 8L2, 9C2, 5P1, 10L2w
	EMED Small Parts	bk3000/5000	13L4w, 18L5, 14L3, 8L2, 18L5s
	EMED Cardiac	bk3000/5000	5P1
EMED OB	bk3000/5000	9C2, E10C4, 6C2, 5C1e	
MSK & Nerve UA2453	MSK & Nerve	bkSpecto*	6C2, 9C2, 5C1e, 5P1e, 13L4w, 14L3, 14L3e, 18L5, X18L5s, 8L2, 18L5s
Women's Health UA2455	OB, GYN	bkSpecto*	6C2, 9C2, 5C1e, E13C2, E10C4
Vascular UA2456	Vascular	bkSpecto*	9C2, 13L4w, 14L3, 14L3e, 18L5, X18L5s, 8L2, 18L5s

	Prostate Triplane E14C4t (9018)	High Resolution Endfire E13C2 (9029)	Endocavity Biplane E14CL4b (9048)	High Resolution Biplane E11C3b (9008)	Endocavity E10C4 (9019)	Endocavity 3D X14L4 (9038)	Anorectal 3D 20R3 (9052)	Linear 14L3 (9051)	Linear 8L2 (9032)	Wide Linear 10L2w (9022)
										
MAIN USES (System Dependent)	Transrectal, Transvaginal	Transrectal, Transvaginal	Transrectal, Transvaginal	Transrectal	Transvaginal, Transrectal, Fetal incl. Obstetrics	Transrectal, Transvaginal	Transrectal, Transvaginal	Small Parts, MSK, PV, Pediatrics	PV, Pediatrics, Small Parts, MSK	PV
Footprint/Tip diameter (mm)	20	20	20	20	22 x 18	20	17	45 x 14	45 x 14	57 x 10
Puncture	•	•	•	•	•			•	•	•
Contrast Imaging	•	•	•			•			•	
Smart Button™	•	•	•	•		•	•	•	•	•
Elastography	•	•	•	•	•	•		•	•	•
VFI								•	•	•
Sterilizable**	•	•	•	•	•	• ^{a)}	• ^{a)}	•	•	•










OPERATES WITH FOLLOWING SYSTEMS

bk5000	•	•	•		•	•	•	•	•	•
bk3000	•	•	•		•	•	•	•	•	•
bkSpecto*	•	•	•	•	•			•	•	

	High Frequency Linear 18L5 (9070)	Small High Frequency Linear 18L5s (9081)	Linear 14L3e (9086)	Wide Linear 13L4w (9011)	High Frequency Curved 9C2 (9002)	Curved 5C1e (9085)	Curved 6C2 (9040)	Small Curved 6C2s (9023)	Small Footprint Cardiac 5P1 (9077)	Small Footprint Cardiac 5P1e (9087)
										
MAIN USES (System Dependent)	Small Parts, MSK, PV	MSK, Small Parts, PV, Pediatrics	Small Parts, MSK, Pediatrics, PV	Small Parts, MSK, PV, Pediatrics	Abdominal, MSK, Fetal incl. Obstetrics, PV, Pediatrics	Abdominal, Obstetrics, Fetal, MSK	Abdominal, MSK, Fetal incl. Obstetrics	Abdominal, Pediatrics, Fetal incl. Obstetrics	Cardiac ^{b)} , Abdominal, Transcranial	Cardiac ^{b)} , Abdominal, Pediatric
Footprint/Tip diameter (mm)	48 x 13	24 x 11	45 x 14	57 x 10	54 x 10	68 x 19	69 x 19	33 x 12	26 x 20	26 x 20
Puncture	•		•	•	•	•	•	•		
Contrast Imaging							•	•		
Smart Button™	•	•	•	•	•	•	•	•	•	•
Elastography	•	•	•	•	•					
VFI			•							
Sterilizable**	•			•	•		•	•	•	

OPERATES WITH FOLLOWING SYSTEMS

bk5000	•	•		•	•	•	•	•	•	
bk3000	•	•		•	•	•	•	•	•	
bkSpecto*	•	•	•	•	•	•	•	•		•

	4-Way Laparoscopic I12C4f (9066)	Intraoperative I14C5I (9015)	T-shaped Intraoperative I14C5T (9016)	Intraoperative Biplane I12C5b (9024)	Robotic Drop-In X12C4 (9026)	Craniotomy N13C5 (9062)	Burr-Hole N11C5s (9063)	"Hockey-Stick" X18L5s (9009)
								
MAIN USES (System Dependent)	Intraoperative	Intraoperative, Pediatrics	Intraoperative, Abdominal	Intraoperative, MSK, PV	Intraoperative	Intraoperative, Neurosurgery, Neonatal, Cephalic, Pediatrics	Intraoperative, Neurosurgery	Intraoperative, PV, Small Parts, MSK, Brain/spinal cord ^{d)}
Footprint/Tip diameter (mm)	33 x 9	60 x 14	60 x 14	40 x 30	33 x 9	29 x 10	10 x 9	31 x 9
Puncture	•	•		•		•	•	•
Contrast Imaging	•	•	•	•	•			
Smart Button™						•	•	•
Elastography ^{d)}	•	•	•	•	•			•
VFI								
Sterilizable**	•	•	•	•	•	•	•	•

OPERATES WITH FOLLOWING SYSTEMS

bk5000	•	•	•	•	•	•	•	•
bk3000						•		•
bkSpecto*								•

a) Sterilization with STERIS System 1 models only.
b) Not for direct use on the heart.
c) Use of the N13C5 and X18L5s for elastography has not been CE-marked.

d) Use of the 'hockey stick' transducer for intraoperative (neuro) has not been CE-marked or licensed by Health Canada

MSK = Musculoskeletal PV = Peripheral Vascular
* bkSpecto has not been licensed by Health Canada
** For complete information about sterilization compatibility and validity, please see Care & Cleaning.