

# DESIGNED

bk5000

# FOR SURGERY

- Superb images instantly - no adjustment needed
- Small system footprint - fits easily in the OR
- Wide range of specialized and sterilizable transducers
- Sterilizable remote control for convenient control of the system in the sterile field
- Rapid boot time and up to two hours of plug-free imaging
- Advanced modes: elastography, contrast imaging, 3D, vector flow imaging (VFI)



## HPB, INTRAOPERATIVE

Thrombus in portal vein  
I12C5b transducer

Elastography of metastatic lesion in liver  
I14C5T transducer

I14C5I Intraoperative  
I14C5T T-Shaped Intraoperative  
I12C5b Intraoperative Biplane  
I12C4f 4-Way Laparoscopic

## NEUROSURGERY, INTRAOPERATIVE

Highly-detailed brain tissue  
N13C5 craniotomy transducer

High-grade tumor  
N13C5 craniotomy transducer

N13C5 High Frequency Curved  
N11C5s Burr-Hole  
X18L5s "Hockey-Stick"

## ROBOTIC-ASSISTED

Renal tumor  
X12C4 Drop-In transducer

X12C4 Robotic Drop-In

## ENDOCAVITY

Expanded sector imaging of prostate gland  
E14C4t transducer

bkFusion Biopsy of region of interest  
E14C4t transducer

E14C4t Prostate Triplane  
E13C2 High Resolution Endfire  
E14C4 Endfire Curved  
E14CL4b Endocavity Biplane  
E10C4 Endocavity

## COLORECTAL, PELVIC FLOOR

T3 rectal mass  
20R3 transducer

20R3 Anorectal 3D  
9038 3D Endocavity

## MINIMALLY INVASIVE

Color Doppler of liver hemangioma  
I12C4f transducer

I12C4f 4-Way Laparoscopic

## PHASED

5P1 Small Footprint Cardiac

## CURVED

6C2 Curved  
5C1e Curved  
6C2s Small Curved  
9C2 High Frequency Curved

## LINEAR

13L4w Wide Linear  
14L3 Linear  
8L2 Linear  
18L5 High Frequency Linear  
18L5s Small High Frequency Linear  
10L2w Wide Linear

N11C5s, 5C1e, E13C2 and 18L5s transducers for bk5000 have not been licensed by Health Canada  
N13C5 for bk5000 is licensed by Health Canada only for surface applications  
5P1 transducer is not for direct use on the heart  
X18L5s for neurosurgery is for future release

BK Ultrasound  
8 Centennial Drive  
Peabody,  
MA 01960  
USA  
T +1 978 326 1300  
bkultrasound.com

USA  
Sales & Service  
BK Ultrasound  
8 Centennial Drive  
Peabody, MA 01960  
USA  
T +1 978 326 1300  
F +1 978 326 1399  
bkultrasound.com

Europe and Rest of World  
Sales, Service & Design Center  
BK Ultrasound  
Mileparken 34,  
2730 Herlev, Denmark  
T +45 4452 8100  
F +45 4452 8199  
bkultrasound.com

Asia  
Sales & Service  
Analogic Medical Equipment  
(Shanghai) Co., Ltd.  
1377, Lan Dian Road  
Pu Dong New District  
Shanghai  
China 201132  
T +86 21 2089 0333  
bkultrasound.com

**bk**  
ultrasound  
powered by analogic